



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

## BILIARY CATHETER



**ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET  
SOFT TIP / RADIOPAQUE**

*INTENDED FOR REMOVAL OF GALL STONES FROM THE BILIARY DUCT*

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
3F	GREEN	3F BIL 23	0.10 CC	5 MM
3F	GREEN	3F BIL 40	0.10 CC	5 MM
4F	RED	4F BIL 23	0.50 CC	9 MM
4F	RED	4F BIL 40	0.50 CC	9 MM
5F	WHITE	5F BIL 23	0.75 CC	11 MM
5F	WHITE	5F BIL 40	0.75 CC	11 MM
6F	BLUE	6F BIL 23	1.25 CC	13 MM
6F	BLUE	6F BIL 40	1.25 CC	13 MM
7F	YELLOW	7F BIL 23	2.00 CC	14 MM
7F	YELLOW	7F BIL 40	2.00 CC	14 MM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: [lmi@earthlink.net](mailto:lmi@earthlink.net)

[www.lucasmedicalinc.com](http://www.lucasmedicalinc.com)

ISO 13485:2016  
REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

GPM S.L.  
GERIFALTE 117  
28232 MADRID, SPAIN  
TEL: +34 91 630 8273  
FAX: +34 91 630 8274



0  
4  
5  
9