




ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE
CAROTID ENDARTERECTOMY SHUNT CATHETER

INTENDED FOR PROVIDING TEMPORARY BLOOD FLOW TO CEREBRAL AREA DURING A CAROTID ENDARTERECTOMY PROCEDURE



- ❖ ALL SILICONE
- ❖ NON-FRAGMENTING BALLOON
- ❖ INTEGRATED FLEXIBLE STYLET
- ❖ SOFT TIP
- ❖ RADIOPAQUE

AVAILABLE SIZES:

MODEL: C B S

FRENCH SIZE	PRODUCT DESIGNATION	RECOMMENDED LIQUID VOLUME (CC)		INFLATED DIAMETER (MM)		LENGTH (CM)
		INTERNAL BALLOON	COMMON BALLOON	INTERNAL BALLOON	COMMON BALLOON	
8F	8F CBS100-10-1	3/4	1 ¼	11	12	15 CM
	8F CBS400-40-1					30 CM
9F	9F CBS100-10-1	3/4	1 ¼	12	13	15 CM
	9F CBS400-40-1					30 CM
10F	10F CBS100-10-1	3/4	1 ¼	12	13.5	15 CM
	10F CBS400-40-1					30 CM
11F	11F CBS100-10-1	3/4	1 ¼	12.5	13.5	15 CM
	11F CBS400-40-1					30 CM
12F	12F CBS100-10-1	3/4	1 ¼	12.5	14	15 CM
	12F CBS400-40-1					30 CM
8F	8F CBS100-10-2	3/4	1 ¼	11	12	15 CM
	8F CBS400-40-2					30 CM
9F	9F CBS100-10-2	3/4	1 ¼	12	13	15 CM
	9F CBS400-40-2					30 CM
10F	10F CBS100-10-2	3/4	1 ¼	12	13.5	15 CM
	10F CBS400-40-2					30 CM
11F	11F CBS100-10-2	3/4	1 ¼	12.5	13.5	15 CM
	11F CBS400-40-2					30 CM
12F	12F CBS100-10-2	3/4	1 ¼	12.5	14	15 CM
	12F CBS400-40-2					30 CM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: lmi@earthlink.net

www.lucasmedicalinc.com

ISO 13485:2016

REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

GPM S.L.
 GERIFALTE 117
 28232 MADRID, SPAIN
 TEL: +34 91 630 8273
 FAX: +34 91 630 8274



0
4
5
9