



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

## ARTERIAL EMBOLLECTOMY CATHETER



**ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET  
SOFT TIP / RADIOPAQUE**

*INTENDED FOR INTRA-OPERATIVE REMOVAL OF ARTERIAL EMBOLI FROM ARTERIES*

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
2F	PURPLE	2F EMB 40	0.05 CC	4 MM
2F	PURPLE	2F EMB 60	0.05 CC	4 MM
3F	GREEN	3F EMB 40	0.10 CC	5 MM
3F	GREEN	3F EMB 80	0.10 CC	5 MM
4F	RED	4F EMB 40	0.50 CC	9 MM
4F	RED	4F EMB 80	0.50 CC	9 MM
4F	RED	4F EMB 100	0.50 CC	9 MM
5F	WHITE	5F EMB 40	0.75 CC	11 MM
5F	WHITE	5F EMB 80	0.75 CC	11 MM
6F	BLUE	6F EMB 40	1.25 CC	13 MM
6F	BLUE	6F EMB 80	1.25 CC	13 MM
7F	YELLOW	7F EMB 40	2.00 CC	14 MM
7F	YELLOW	7F EMB 80	2.00 CC	14 MM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: [lmi@earthlink.net](mailto:lmi@earthlink.net)

[www.lucasmedicalinc.com](http://www.lucasmedicalinc.com)

ISO 13485:2016

REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

GPM S.L.  
GERIFALTE 117  
28232 MADRID, SPAIN  
TEL: +34 91 630 8273  
FAX: +34 91 630 8274



0  
4  
5  
9