



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

THROMBECTOMY CATHETER



ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET
SOFT TIP / RADIOPAQUE

INTENDED FOR INTRA-OPERATIVE REMOVAL OF VENOUS THROMBOLI FROM VEINS

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
4F	RED	4F THB 80	0.50 CC	9 MM
5F	WHITE	5F THB 80	0.75 CC	11 MM
6F	BLUE	6F THB 80	1.25 CC	13 MM
7F	YELLOW	7F THB 80	2.00 CC	14 MM
8F	BROWN	8F THB 80	3.00 CC	16 MM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: lmi@earthlink.net

www.lucasmedicalinc.com

ISO 13485:2016
REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

GPM S.L.
GERIFALTE 117
28232 MADRID, SPAIN
TEL: +34 91 630 8273
FAX: +34 91 630 8274

