



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

BILIARY CATHETER



ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET
SOFT TIP / RADIOPAQUE

INTENDED FOR REMOVAL OF GALL STONES FROM THE BILIARY DUCT

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
3F	GREEN	3F BIL 23	0.10 CC	5 MM
3F	GREEN	3F BIL 40	0.10 CC	5 MM
4F	RED	4F BIL 23	0.50 CC	9 MM
4F	RED	4F BIL 40	0.50 CC	9 MM
5F	WHITE	5F BIL 23	0.75 CC	11 MM
5F	WHITE	5F BIL 40	0.75 CC	11 MM
6F	BLUE	6F BIL 23	1.25 CC	13 MM
6F	BLUE	6F BIL 40	1.25 CC	13 MM
7F	YELLOW	7F BIL 23	2.00 CC	14 MM
7F	YELLOW	7F BIL 40	2.00 CC	14 MM

1554 NORTH CASE STREET, ORANGE, CA 92867 PHONE: (714) 938-0233

E-MAIL: lmi@earthlink.net www.lucasmedicalinc.com

ISO 13485:2016
REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332