



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

ARTERIAL EMBOLICATHETER



ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET
SOFT TIP / RADIOPAQUE

INTENDED FOR INTRA-OPERATIVE REMOVAL OF ARTERIAL EMBOLI FROM ARTERIES

AVAILABLE SIZES:

| FRENCH SIZE | COLOR CODE | PRODUCT DESIGNATION | MAXIMUM LIQUID INFLATION VOLUME (MLV) | INFLATED BALLOON DIAMETER @ MLV |
|-------------|------------|---------------------|---------------------------------------|---------------------------------|
| 2F | PURPLE | 2F EMB 40 | 0.05 CC | 4 MM |
| 2F | PURPLE | 2F EMB 60 | 0.05 CC | 4 MM |
| 3F | GREEN | 3F EMB 40 | 0.10 CC | 5 MM |
| 3F | GREEN | 3F EMB 80 | 0.10 CC | 5 MM |
| 4F | RED | 4F EMB 40 | 0.50 CC | 9 MM |
| 4F | RED | 4F EMB 80 | 0.50 CC | 9 MM |
| 4F | RED | 4F EMB 100 | 0.50 CC | 9 MM |
| 5F | WHITE | 5F EMB 40 | 0.75 CC | 11 MM |
| 5F | WHITE | 5F EMB 80 | 0.75 CC | 11 MM |
| 6F | BLUE | 6F EMB 40 | 1.25 CC | 13 MM |
| 6F | BLUE | 6F EMB 80 | 1.25 CC | 13 MM |
| 7F | YELLOW | 7F EMB 40 | 2.00 CC | 14 MM |
| 7F | YELLOW | 7F EMB 80 | 2.00 CC | 14 MM |

1554 NORTH CASE STREET, ORANGE, CA 92867 PHONE: (714) 938-0233

E-MAIL: lmi@earthlink.net www.lucasmedicalinc.com

ISO 13485:2016
REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

GPM S.L.
GERIFALTE 117
28232 MADRID, SPAIN
TEL: +34 91 630 8273
FAX: +34 91 630 8274



0
4
5
9