



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

THROMBECTOMY CATHETER



**ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET
SOFT TIP / RADIOPAQUE**

INTENDED FOR INTRA-OPERATIVE REMOVAL OF VENOUS THROMBOLI FROM VEINS

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
4F	RED	4F THB 80	0.50 CC	9 MM
5F	WHITE	5F THB 80	0.75 CC	11 MM
6F	BLUE	6F THB 80	1.25 CC	13 MM
7F	YELLOW	7F THB 80	2.00 CC	14 MM
8F	BROWN	8F THB 80	3.00 CC	16 MM

1554 NORTH CASE STREET, ORANGE, CA 92867 PHONE: (714) 938-0233
E-MAIL: lmi@earthlink.net www.lucasmedicalinc.com

ISO 13485:2016
REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332